HAMAMATSU PHOTOMULTIPLIER TUBES R4220 R4220P (For Photon Counting)

Very High Cathode Sensitivity with Low Noise Photocathode

FEATURES

Spectral Response	. 185 nm to 710 nm
High Cathode Sensitivity	
Luminous	100 μA/Im
Radiant at 410 nm	70 mA/W
High Anode Sensitivity (at 1000 V)	
Luminous	1200 A/Im
Radiant at 410 nm	8.4 × 10 ⁵ A/W
Low Dark Current	0.2 nA
●Low Dark Counts (R4220P)	10 s ⁻¹

APPLICATIONS

Fluorescence Spectrometer
Chemiluminescence Detection
Raman Spectroscopy
Low Light Level Ditection




SPECIFICATIONS

GENERAL

Parameter	Description / Value	Unit	
Spectral Response	185 to 710	nm	
Wavelength of Maximum Response	410	nm	
Photocathode			
Material	Low noise bialkali	—	
Minimum Effective Area	8×24	mm	
Window Material	UV glass	-	
Dynode			
Secondary Emitting Surface	Low noise bialkali	—	
Structure	Circular-cage	—	
Number of Stages	9	—	
Direct Interelectrode Capacitances			
Anode to Last Dynode	4	pF	
Anode to All Other Electrodes	6	рF	
Base	11-pin base	_	
	JEDEC No. B11-88		
Weight	45	g	
Operating Ambient Temperature	-30 to +50	°C	
Storage Temperature	-30 to +50	°C	
Suitable Socket	E678–11A (Sold Separately)	_	
Cuitable Cooket Assembly	E717–63 (Sold Separately)		
Suitable Socket Assembly	E717-74 (Sold Separately)		

Figure 1: Typical Spectral Response



Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office. Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2010 Hamamatsu Photonics K.K.

PHOTOMULTIPLIER TUBES R4220, R4220P (For Photon Counting)

MAXIMUM RATINGS (Absolute Maximum Values)

Parameter	Value	Unit
Supply Voltage		
Between Anode and Cathode (DC)	1250	V
Between Anode and Last Dynode (DC)	250	V
Average Anode Current [®]	0.1	mA

CHARACTERISTICS (at 25 °C)

R4220 R4220P for General Purpose for Photon Counting Unit **Parameter** Min. Typ. Max. Min. Typ. Max. Cathode Sensitivity 23 23 % Quantum Efficiency (at peak wavelength) Luminous [®] 80 100 80 100 μA/lm Radiant (at peak wavelength) 70 _ 70 mA/W Blue Sensitivity Index © 8 8 ____ Anode Sensitivity Luminous D 1000 1200 1000 1200 A/Im Radiant at 400 nm 8.4×10^{5} 8.4×10^{5} A/W Gain ® 1.2×10^{7} 1.2×10^{7} Anode Dark Current © After 30 minutes Storage in the darkness 0.2 0.2 2.0 0.5 nA Anode Dark Counts (E) 10 50 S-1 ENI(Equivalent Noise Input) © 3.30×10^{-17} 3.30×10^{-17} _ _ W Time Response ^(D) Anode Pulse Rise Time ^(II) 2.2 2.2 ns Electron Transit Time J 22 22 ns Transit Time Spread (TTS) ® 1.2 1.2 ns Anode Current Stability © **Current Hysteresis** 0.1 0.1 % Voltage Hysteresis % 1.0 1.0

NOTES

- A: Averaged over any interval of 30 seconds maximum.
- B: The light source is a tungsten filament lamp operated at a distribution temperature of 2856 K. Supply voltage is 150 volts between the cathode and all other electrodes connected together as anode.
- C: The value is cathode output current when a blue filter(Corning CS-5-58 polished to 1/2 stock thickness) is interposed between the light source and the tube under the same condition as Note B.
- D: Measured with the same light source as Note B and with the anode-tocathode supply voltage and voltage distribution ratio shown in Table 1 below.
- E: Measured with the same supply voltage and voltage distribution ratio as Note D after removal of light.
- F: Measured at the plateau voltage.
- G:ENI is an indication of the photon-limited signal-to-noise ratio. It refers to the amount of light in watts to produce a signal-to-noise ratio of unity in the output of a photomultiplier tube.


$$\mathsf{ENI} = \frac{\sqrt{2q} \cdot \mathsf{Idb} \cdot \mathsf{G} \cdot \mathsf{f}}{\mathsf{S}}$$

where $q = Electronic charge (1.60 \times 10^{-19} coulomb).$

- Idb = Anode dark current(after 30 minute storage) in amperes. G = Gain.
 - f = Bandwidth of the system in hertz. 1 hertz is used.
 - S = Anode radiant sensitivity in amperes per watt at the wavelength of peak response.
- H: The rise time is the time for the output pulse to rise from 10% to 90% of the peak amplitude when the entire photocathode is illuminated by a delta function light pulse.
- J: The electron transit time is the interval between the arrival of delta function light pulse at the entrance window of the tube and the time when the anode output reaches the peak amplitude. In measurement, the whole photocathode is illuminated.

K: Also called transit time jitter. This is the fluctuation in electron transit time between individual pulses in the signal photoelectron mode, and may be defined as the FWHM of the frequency distribution of electron transit times.

L: Hysteresis is temporary instability in anode current after light and voltage are applied.



(1)Current Hysteresis

The tube is operated at 750 volts with an anode current of 1 micro-ampere for 5 minutes. The light is then removed from the tube for a minute. The tube is then re-illuminated by the previous light level for a minute to measure the variation.

(2)Voltage Hysteresis

The tube is operated at 300 volts with an anode current of 0.1 micro-ampere for 5 minutes. The light is then removed from the tube and the supply voltage is quickly increased to 800 volts. After a minute, the supply voltage is then reduced to the previous value and the tube is re-illuminated for a minute to measure the variation.

Table 1:Voltage Distribution Ratio

Electrodes	ł	<	Dy1	Dy	'2 C)уЗ	Dy	/4	Dy5	Dy	/6	Dy	7	Dy	/8	Dy	/9		Р
Distribution Ratio		1		1	1		1	1		1		1		1		1		1	

Supply Voltage : 1000 V (DC)

K : Cathode, Dy : Dynode, P: Anode

HAMAMATSU

Figure 2: Typical Gain and Anode Dark Current




Figure 4: Typical ENI vs. Wavelength




Figure 3: Typical Time Response





Figure 5: Typical EADCI (Equivalent Anode Dark Current Input) vs. Supply Voltage



Data shown here, which is given from a relation among supply voltage, anode sensitivity and dark current, serves as a good reference in order to determine the most suitable supply voltage or its range.

Figure 6: Typical Plateau Data for R4220P







Figure 7: Typical Temperature Characteristics of Dark Count for R4220P



PHOTOMULTIPLIER TUBES R4220, R4220P (For Photon Counting)


Figure 8: Dimensional Outline and Basing Diagram (Unit: mm)






TACCA0064EA









* Hamamatsu also provides C4900 series compact high voltage power supplies and C6270 series DP type socket assemblies which incorporate a DC to DC converter type high voltage power supply.

Warning–Personal Safety Hazards

Electrical Shock–Operating voltages applied to this device present a shock hazard.

HAMAMATSU

HAMAMATSU PHOTONICS K.K., Electron Tube Division

WEB SITE www.hamamatsu.com

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205 U.S.A.: Hamamatsu Corporation: 360 Foothill Road, P. O. Box 6910, Bridgewater. N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Hersching am Ammersee, Germany, Telephone: (4)9152-375-0, Fax: (4)9152-2658 E-mail: info@hamamatsu.com France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saulie Trapu, Pare du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (34)8152-375-0, Fax: (4)9152-2658 E-mail: info@hamamatsu.de France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saulie Trapu, Pare du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (34)8 65 371 10, Fax: (34)105 5371 10, Fax: (40)1707-326777 E-mail: info@hamamatsu.co.uk North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171-41 SOLNA, Sweden, Telephone: (46)8-509-031-01, E-mail: info@hamamatsu.se Hawamatsu Photonics Italia: S.R.L.: Strada della Moia, 1/E, 2020 Arese, (Milano), Italy, Telephone: (36)02-935 81 731, Fax: (39)02-935 81 741 E-mail: info@hamamatsu.it

TPMS1003E03 APR. 2010. IP

Figure 9: Socket E678-11A (Sold Separately)