

# PMA-20

Photonic multichannel analyzer C10494 series

Continuous spectrum measurement with 100  $\mu$ s temporal resolution



The PMA-20 Photonic multichannel analyzer C10494 series measures spectral change of emission, absorption and reflection with high sensitivity 100  $\mu$ s temporal resolution.

The PMA-20 is capable of measuring a wide spectrum range from 200 nm to 1100 nm with high wavelength resolution within 3 nm. Using an optional fiber for light collection, spectrum can be obtained easily by directing the fiber close to the sample. As the wavelength axis and the spectral response characteristic are calibrated at the factory, spectral measurements can be carried out easily and accurately.

## Features

- Single shot spectrum measurement with 100  $\mu$ s temporal resolution
- High speed integration
- A compact unit with a spectrometer, a photo-detector and a power supply
- Easy measurements by optical fiber
- Factory calibrated spectral response and wavelength axis characteristics


## Applications

- Time resolved spectrum measurement for emission
- Protein-protein interaction analysis with absorption spectrum
- Chemical reaction tracking with a stopped-flow method
- Photo physics and laser spectroscopy with submillisecond temporal resolution

Example of Measurement


1 Absorption spectrum measurement of BTB solution (pH indicator)

● Absorption measurement system




The PMA-20 measures the spectrum change of the weak acid BTB solutions from yellow to blue by adding sodium bicarbonate with 100 μs temporal resolution.

● Absorption spectrum measurement of BTB solution (pH indicator)




2 Spectrum change measurement with stopped-flow method

● Stopped-flow method measurement system



The PMA-20 measures the reproduction process of rhodopsin as opsin and retinal are mixed with the stopped-flow method.


● Change in absorption



\* Data is provided by courtesy of Department of biophysics Division of biological sciences Graduate school of science Kyoto university.


3 Phosphorescence measurement

● Phosphorescence measurement system



The PMA-20 measures the attenuation process of phosphorescence materials excited by YAG laser in single shot. It also measures by single shot to excite the materials with YAG laser (266 nm, 3 mJ).

● Phosphorescence measurement



**Options**



**Sample Holder for transmission and fluorescence measurement A6751**

This is a dedicated holder with an integrated condensing lens for the use with vials.



**Reflection measurement optics A9665**

These are optics making it possible to illuminate the sample at 45° from the light source and measure the reflected light.



**Variable angle reflection measure optics A10687**

These are optics making it possible to change the angle of input and output ports at maximum 65° and measure the reflected light and fluorescence.



**Digital delay generator C13430-01**

This outputs the gate pulse necessary for an external trigger and gate operation.



**C-mount fiber adapter A6399**

This is an adapter for securing the fiber input optics to the C-mount of a micro-scope or the like.



**C-mount adapter for positioning A9607**

In addition to the function of the C-mount fiber adapter, the measurement position can be checked. Measurements in the UV.



**Attenuation fiber adapter A10474-01**

This adaptor is used when the light power is too strong. It can reduce the input light power by using a pin-hole. (fading rate approx. 1/20 to 1/500)



**OBJECTIVE LENS A4869**

Condensing lens for UV. f=50 mm, F3.5 (A6399, A8482 required)




**Integrating sphere A5640**

This is the integrating sphere for getting complete diffuse light. You can get even intensity light without spread of light source or influence of directional characteristics. (A6399 required)




**Halogen lamp L6758-11**

This is a halogen light source with output wavelengths from 400 nm to 1600 nm for excitation and absorption measurements.



**Xe light source High stability 150 W L6759**

This is a high stability xenon light source with output wavelengths from 250 nm to 1600 nm for excitation and absorption measurements.



**Data analyzer Desktop type C10471-12**

Collecting and analyzing data.

OS	Windows 7 (64 bit)
Memory	2 GB or more
Hard disk	80 GB or more
Interface	Camera Link



Minilite-II Surelite-II-10

**Nd : YAG Laser**

Nd: YAG laser is the excitation light source for transient absorption measurements. The system can be enhanced by combining with the lasers from various manufactures.

**Recommended laser**

- From Continuum : Minilite-II  
Surelite-II-10
- From Spectra-Physics, Inc., : INDI-40




INDI-40

\* For details, please consult the Hamamatsu's local office of each country.

**Specifications**

Type number	C10494-01	C10494-02
Photo-detector	BT-CCD linear image sensor	
Wavelength	200 nm to 950 nm	350 nm to 1100 nm
Wavelength resolution (FWHM)	<3 nm	
Exposure time	0.1 ms to 1 s	
Number of photosensitive device channels	2048 ch	
Pixel size	12 μm × 972 μm	
Read-out noise	100 electrons	
Dark current	100 electrons/scan (+25 °C 20 ms)	
AD resolution	12 bit	
Spectrograph	Czerny-turner type	
Spectrograph F number	4	
Fiber	Quartz fiber 1.5 m	
Interface	Camera Link, USB 2.0	
Power supply	AC100 V to AC240 V, 50 Hz/60 Hz	



**Software**


■ Basic software for PMA-20 U6039-08

- Measurement functions ..... Monitoring measurement  
Data measurement
- Time resolved spectrum measurement ... Emission spectra  
Reflectance  
Transmittance
- Data acquisition condition settings ... Exposure time settings  
Memory integration times setting
- Calibration/correction ..... Wavelength axis calibration  
Sensitivity inconsistency calibration  
Dark current correction
- Display functions ..... Spectrum display  
Time profile display
- Wavelength axis display ..... Wavelength, Wave number, Raman shift, Energy (eV)
- Brightness axis display ..... Linear, Logarithm
- Cursor analysis functions ..... Wavelength (wave number, etc.) vs. Intensity  
Peak detection  
FWHM measurement  
Integrated intensity
- Other analytical functions ..... Smoothing  
Differential waveform operation  
Color calculation (XYZ, xy, uv, Lab)
- Trigger functions ..... Internal trigger, Opening trigger,  
External exposure start, Optical trigger

- Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.
  - Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
  - Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.
- Specifications and external appearance are subject to change without notice.

© 2016 Hamamatsu Photonics K.K.

**Dimensional Outlines**



**HAMAMATSU PHOTONICS K.K. www.hamamatsu.com**

**HAMAMATSU PHOTONICS K.K., Systems Division**

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com  
 Germany: Hamamatsu Photonics Deutschland GmbH.: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de  
 France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr  
 United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk  
 North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se  
 Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-935-81-733, Fax: (39)02-935-81-741 E-mail: info@hamamatsu.it  
 China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn  
 Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No.158, Section2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)03-659-0080, Fax: (886)07-811-7238 E-mail: info@tw.hpk.co.jp