

Head-on PMT

Photon Counting Head H9319 Series



The H9319 series photon counting heads are designed to perform photon counting by just connecting to a PC. The H9319 series includes a 25-mm (1") diameter head-on photomultiplier tube, a photon counting circuit, a high-voltage power supply circuit, counter and a microprocessor. Data transfer, measurement time and other necessary adjustments can be controlled by commands from the PC through the RS-232C interface.

Since the H9319 series performs linearity correction by the internal microprocessor, it provides excellent count linearity within a range of $\pm 1\%$ at $20 \times 10^6 \text{ s}^{-1}$.

Product Variations

| Type No. | Spectral Response | Sample Program | Prescaler |
|----------|-------------------|----------------|-----------|
| H9319-01 | 300 nm to 650 nm | yes | 1/4 |
| H9319-11 | | no | 1/4 |
| H9319-02 | 300 nm to 850 nm | yes | 1/4 |
| H9319-12 | | no | 1/4 |

This product can't be used at vacuum environment or reduced pressure environment.

Specifications

(at +25 °C)

| Parameter | | H9319-01 | H9319-11 | H9319-02 | H9319-12 | Unit |
|------------------------------------------|------|---------------------------------------------------------|-------------------|----------|-------------------|--------------------------------------|
| Input Voltage | | +4.75 to +5.25 | | | | V |
| Max. Input Voltage | | +6 | | | | V |
| Max. Input Current | | 60 | | | | mA |
| Effective Area | | $\phi 22$ | | | | mm |
| Peak Sensitivity Wavelength | | 420 | | | | nm |
| Count Sensitivity | Typ. | 300 nm | 2.3×10^5 | | 2.1×10^5 | $\text{s}^{-1} \cdot \text{pW}^{-1}$ |
| | | 400 nm | 4.1×10^5 | | 2.5×10^5 | |
| | | 500 nm | 3.4×10^5 | | 2.0×10^5 | |
| | | 600 nm | 5.7×10^4 | | 1.3×10^5 | |
| | | 700 nm | — | | 7.8×10^4 | |
| Count Linearity *1 | | 20×10^6 | | | | s^{-1} |
| Dark Count *2 | Typ. | 150 | | 10 000 | | s^{-1} |
| | Max. | 300 | | 15 000 | | |
| PMT Operating Voltage Range | | +300 to +1200 | | | | V |
| Integration Time | | 10 to 1000 | | | | ms |
| Settling Time | | 1 *3 | | | | s |
| | | 5 *4 | | | | s |
| Input Signal (External Trigger Input) *5 | | TTL level signal | | | | — |
| Output Signal (User Line Output) *6 | | TTL level signal | | | | — |
| Interface | | RS-232C, 9600 baud, Parity none, 8 data bit, 1 stop bit | | | | — |
| Operating Ambient Temperature *7 | | +5 to +50 | | | | °C |
| Storage Temperature *7 | | -20 to +50 | | | | °C |
| Weight *8 | | 280 | | | | g |

*1: Random pulse, within $\pm 1\%$ count loss (by count linearity compensation)

*2: After 30 minutes storage in darkness

*3: The time required for the output to reach a stable level following a change in the control voltage from 500 V to 1000 V in darkness

*4: The time required for the output to reach a stable level following a change in the control voltage from 1000 V to 500 V in darkness

*5: Selectable ECR mode or LCR mode (refer to command list)


*6: Controllable by RS-232C command

*7: No condensation


*8: Main body: Approx. 180 g

Photon Counting Head with CPU+Interface

Characteristic (Count sensitivity)



Block Diagram



Command List

| Action | Command *9 | Explanation | Argument | Response *10 |
|---------------------------------|-------------------|------------------------------------------------------------------------------------------------------------------|-------------|----------------|
| Set the Integration Time | P#C _R | Set the number of 10 msec intervals to sum. Same as integration time for 1 reading. | 1 to 100 | VA, BC, BA |
| Set the Sequence of Readings | R#C _R | Set sequence of readings, where each reading uses the integration time set with the P command. | 1 to 255 | VA, BC, BA |
| Set the PMT Input Voltage | V##C _R | Change the high voltage applied to the tube. | 300 to 1200 | VA, BC, BA |
| | D C _R | Re-set the default high voltage to the tube | — | VA, BC |
| Set the Output of the User Line | O#C _R | Set the output of the user digital output line. 0: Low level output 1: High level output | 0 or 1 | VA, BC, BA |
| | S C _R | Start the reading sequence | — | 4 byte/reading |
| Start the Reading Sequence | C C _R | Start a continuous reading process. Will continue indefinitely until a STOP character (C _R) is sent. | — | 4 byte/reading |
| | E C _R | Start a reading sequence for each positive-edge TTL transition applied to the purple stripe user line. | — | 4 byte/reading |
| | L C _R | Start a reading sequence for each positive-level TTL transition applied to the purple stripe user line. | — | 4 byte/reading |

*9: The Response acknowledgment is returned having two bytes.

*10: VA: valid command, BC: bad command, BA: bad argument

Dimensional Outlines (Unit: mm)

