

Improved detection limit for near infrared (to 900 nm)

## FEATURES

- Low dark current (1/50 of that at room temperature)
- Fast cooling (3 minutes) by thermoelectric cooler directly coupled to the photocathode
- Wide dynamic range
- Free of condensation

## APPLICATIONS

- NO<sub>x</sub> Gas Detection
- Fluorescence Detection (LIF, Fluorescence Spectrophotometer)
- Chemiluminescence Detection
- NIR Spectroscopy



## COOLING SPECIFICATIONS


Parameter	Description / Value	Unit
Cooling Method	Thermoelectric cooling	—
Max. Cooling Temperature ( $\Delta T$ ) *1 *2	20	°C
Cooling Time *1	Approx. 3	min

\*1: Input current to thermoelectric cooler = 2.1 A


\*2: Photocathode temperature difference from ambient

## COOLING CHARACTERISTICS

### ● S/N Ratio During Cooling



### ● Anode dark current after cooling has started



# PHOTOSENSOR MODULE WITH THERMOELECTRIC COOLER H7844

## SPECIFICATIONS

(at 25 °C)

Parameter		Value	Unit
Input Voltage		+11.5 to +15.5	V
Max. Input Voltage for Main Unit		+18	V
Max. Input Current for Main Unit		40	mA
Max. Input Voltage for Thermoelectric Cooler		3.75	V
Max. Input Current for Thermoelectric Cooler		3.9	A
Max. Output Signal Current		58	μA
Max. Control Voltage		+1.2 (Input impedance: 100 kΩ)	V
Recommended Control Voltage Adjustment Range		+0.3 to +1.1	V
Effective Area		10 × 14	mm
Spectral Response Range		185 to 900	nm
Peak Sensitivity Wavelength		400	nm
Cathode	Luminous Sensitivity	Min. 140 Typ. 300	μA/lm
	Blue Sensitivity Index (CS 5-58)	Typ. 9.0	
	Red / White Ratio (R-68)	Typ. 0.3	—
	Radiant Sensitivity *3	Typ. 76	mA/W
Anode	Luminous Sensitivity *4	Min. 400 Typ. 3000	A/lm
	Radiant Sensitivity *3 *4	Typ. $7.6 \times 10^5$	
	Gain *4	Typ. $1.0 \times 10^7$	—
	Dark Current *4 *5	Typ. 0.1 Max. 1.0	nA
Equivalent Noise Input (ENI) *3 *4 *5		$2.4 \times 10^{-17}$	W
Rise Time *4		2.2	ns
Settling Time *6		0.2	s
Operating Ambient Temperature *7		+5 to +40	°C
Storage Temperature *7		-20 to +50	°C
Weight		296	g

\*3: At peak sensitivity wavelength

\*4: Control voltage +1.0 V (PMT supply voltage -1000 V), with cooler operated


\*5: After 30 minutes storage in darkness

\*6: The time required for the output to reach a stable level following a change in the control voltage from +1.0 V to +0.5 V.


\*7: No condensation

## CHARACTERISTICS


### ●Spectral Response




### ●Equivalent Noise Input (ENI)



●Gain vs. Control Voltage




●Plateau Characteristic




**SENSITIVITY ADJUSTMENT AND COOLING OPERATION**

VOLTAGE PROGRAMMING



- Adjust the control voltage when adjusting the anode sensitivity of the PMT.
- Electrically isolate the reference voltage output. (This output is not used.)
- Low voltage input can also be used to supply the power to the cooling fan. In this case, the low voltage input must be set to +12 V.
- Always run the cooling fan while the thermoelectric cooler is operating.

RESISTANCE PROGRAMMING




- Monitor the control voltage when adjusting the anode sensitivity of the PMT with a trimmer potentiometer.
- A C-mount adaptor is available for H7844. (Sold separately)


# PHOTOSENSOR MODULE WITH THERMOELECTRIC COOLER H7844

## DIMENSIONAL OUTLINE (Unit: mm)

\*Dimensional tolerance is  $\pm 0.5$  mm unless otherwise specified.



## ACCESSORIES (Supplied)



TPMOA0036EB

**HAMAMATSU PHOTONICS K.K.** [www.hamamatsu.com](http://www.hamamatsu.com)

**HAMAMATSU PHOTONICS K.K., Electron Tube Division**

**314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205**

**U.S.A.:** Hamamatsu Corporation, 360 Foothill Road, P. O. Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: [usa@hamamatsu.com](mailto:usa@hamamatsu.com)

**Germany:** Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: [info@hamamatsu.de](mailto:info@hamamatsu.de)

**France:** Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: [infos@hamamatsu.fr](mailto:infos@hamamatsu.fr)

**United Kingdom:** Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road Welwyn Garden City Hertfordshire AL7 1BW, United Kingdom, Telephone: 44-(0)1707-294888, Fax: 44(0)1707-325777 E-mail: [info@hamamatsu.co.uk](mailto:info@hamamatsu.co.uk)

**North Europe:** Hamamatsu Photonics Norden AB: Thorshamnsgatan 35 SE-164 40 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: [info@hamamatsu.se](mailto:info@hamamatsu.se)

**Italy:** Hamamatsu Photonics Italia: S.R.L.: Strada della Moia, 1/E, 20020 Arese, (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741 E-mail: [info@hamamatsu.it](mailto:info@hamamatsu.it)

**China:** Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Road North, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: [hpc@hamamatsu.com.cn](mailto:hpc@hamamatsu.com.cn)

TPMO1032E04  
SEPT. 2013 IP